

# Add Taq DNA Polymerase

Research Use Only

## Product Code

17001X

## Component

1. Add Taq DNA Polymerase (5.0 U/  $\mu$ l) 1,000 units, 200  $\mu$ l
2. 10x Reaction Buffer (without Mg) 1.0 ml x 2 tubes
3. 25mM MgCl<sub>2</sub> 1.0 ml x 2 tubes

## Storage Condition

Store at -20°C

## Description

Add Taq DNA Polymerase is a highly thermostable recombinant DNA polymerase derived from the thermophile, *Thermus aquaticus*.

The molecular weight of the recombinant protein is 94kD. Add Taq DNA polymerase catalyzes the 5'→3' synthesis of DNA but has no detectable 3'→5' proofreading exonuclease activity, and possesses low 5'→3' exonuclease activity, which results in a 3'-dA overhang on the PCR product.

## Storage Buffer

20mM Tris-HCl (pH 8.0), 100mM KCl, 3mM MgCl<sub>2</sub>, 1mM DTT, 0.1% Nonidet P-40, 0.1% Tween® 20 and 50% (v/v) glycerol

## 10X Reaction Buffer

100mM Tris-HCl (pH8.8), 500mM KCl, 1% Triton® X-100 and 20mM MgCl<sub>2</sub>

## Storage and Stability

Add Taq DNA Polymerase is stable for 2 years when stored in a constant temperature freezer at less than -20°C.

## Nucleic Acid Amplification Protocol

1. Add the following components to a thin-walled PCR tube:

Nuclease-Free Water	x $\mu$ l
10x Reaction Buffer	2 $\mu$ l
10mM dNTP Mixture	2 $\mu$ l
Forward primer (10 $\mu$ M)	0.25~2 $\mu$ l
Reverse primer (10 $\mu$ M)	0.25~2 $\mu$ l
DNA template	x $\mu$ l
Add Taq DNA Polymerase (5 U/ $\mu$ l)	0.2 $\mu$ l
Total reaction volume	20 $\mu$ l

\* Recommendation for template DNA concentration in a 20  $\mu$ l reaction volume

- 1) Human genomic DNA: 0.1 ng ~ 1  $\mu$ g
- 2) Bacterial genomic DNA: 0.1 ng ~ 100 ng
- 3) Plasmid DNA: 0.01 ng ~ 5 ng

2. PCR cycling

Initial denaturation	95°C, 5 min
	95°C, 15 – 30 sec
PCR cycling (25 – 40 cycles)	55 - 65°C, 15 – 30 sec
	72°C, 30 sec per kb of product length
Final extension	72°C, 5 min
Hold	12°C, $\infty$