

AddScript RT Master (2x Conc.)

Research Use Only

Product Code

22101

Component

1. AddScript RT Master (2x conc.) 1.0 ml

Storage Condition

Store at -20°C

Description

AddScript RT Master provides sensitive and easy-to-use components which contain all the reagents for first strand cDNA synthesis.

Especially, thermostable MMLV RTase (RNase H-), reaction buffer, dNTP mixture, and enzyme stabilizer are included in 2x AddScript RT Master.

Usage Information

- The reaction temperature for cDNA synthesis is 50°C.
- The reaction time for cDNA synthesis is 60 min.

Quality Control

The performance of AddScript RT Master (2x) is tested in an RT reaction using human total RNA with oligo dT₂₀ and random hexamer each. The sensitivity of the kit is verified by the detection of GAPDH and Actin transcript in 10 pg total RNA after 30 cycles.

Storage and Stability

AddScript RT Master (2x) is stable for 2 years when stored in a constant temperature freezer at less than -20°C.

Made in KOREA

This product was manufactured through ISO 9001 & 13485 system.

Reaction Assembly

1. Add the following components to a thin-walled PCR tube:

Nuclease-free D.W	x μ l
2x AddScript RT Master	10 μ l
10mM dNTP Mixture (Not provided)	2.0 μ l
50~100 pmoles/ μ l oligo dT ₂₀ (random hexamer) or Gene specific primer (10~20 pmoles/ μ l)	1.0 μ l
RNA template	x μ l
RNase Inhibitor (Optional)	x μ l
Total reaction volume	20 μ l

* Recommendation for template RNA concentration in a 20 μ l reaction volume

1) total RNA: 100 fg ~ 1 μ g

2) mRNA: 10 fg ~ 1 μ g

2. Temperature cycling Protocol

Priming	25°C, 10 min
Reverse transcription	50°C, 60 min
RT inactivation	80°C, 5 min
Hold	12°C, ∞

Recommendation for downstream PCR

For downstream PCR amplification, the volume of cDNA product should not exceed 1/5 of the PCR reaction volume, typically 1 ~ 4 μ l in 20 μ l PCR reaction.