

AddScript RT-PCR Master (2x Conc.)

Research Use Only

Product Code

24701

Component

1. AddScript RT-PCR Master 1ml

Storage Condition

Store at -20°C

Description

AddScript RT-PCR Master provides sensitive and easy-to-use components which contain all the reagents for first strand cDNA synthesis and PCR reaction in one-tube reaction continuously.

Especially, thermostable MMLV RTase (RNase H-), hot-start Taq DNA Polymerase and RNase Inhibitor, loading dye and sediment for electrophoresis and dNTP mixture are included in 2x Master Mix.

Quality Control

The performance of AddScript RT-PCR Master is tested in a RT and PCR one-tube reaction using human total RNA with specific primers. The sensitivity of the kit is verified by the detection of GAPDH and Actin transcript in 10 pg total RNA after 30 cycles.

Storage and Stability

AddScript RT-PCR Master is stable for 2 year when stored in a constant temperature freezer at less than -20°C.

Made in KOREA

This product was manufactured through ISO 9001 & 13485 system.

Nucleic Acid Amplification Protocol

1. Add the following components to a thin-walled PCR tube:

Nuclease-free D.W	x μ l
2x AddScript RT-PCR Master	10.0 μ l
Forward primer (10 μ M)	0.25~2.0 μ l
Reverse primer (10 μ M)	0.25~2.0 μ l
RNA template	x μ l
Total reaction volume	20 μ l

* Recommendation for template RNA concentration in a 20 μ l reaction volume

1) total RNA: 100 fg ~ 1 μ g

2) mRNA: 10 fg ~ 1 μ g

2. Temperature cycling Protocol

cDNA synthesis	50°C, 30 min
Initial Denaturation	95°C, 10 min
PCR Cycling (30– 40 cycles)	95°C, 15 – 30 sec
	55 - 65°C, 15 – 30 sec
	72°C, 1 min
Final Extension	72°C, 5 min
Hold	12°C, ∞