

AGAROSE LE (Biotechnology Grade II)

Research Use Only

Product Code

A1710 (100 gram)

A1750 (500 gram)

Analytical Specifications

Gelling temperature (1% Gel)	$88^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$
Melting temperature (1% Gel)	$36^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$
Gel Strength (1% Gel)	$\geq 1200\text{g/cm}^2$
EEO (m _r)	≤ 0.13
Sulfate	$\leq 0.15\%$

Storage Condition

Store at room temperature

Storage and Stability

Agroase LE is stable for 3 years when stored in a constant temperature without moisture condition.

Made in KOREA

This product was manufactured through ISO 9001 & 13485 system.